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Vol. 32, No. 8

CHROMagar Candida, a New Differential Isolation Medium for Presumptive Identification of Clinically Important *Candida* Species

FRANK C. ODDS* AND RIA BERNAERTS

Department of Bacteriology and Mycology, Janssen Research Foundation, B-2340 Beerse, Belgium

Received 7 February 1994/Returned for modification 31 March 1994/Accepted 21 April 1994


VOL. 32, 1994

NEW DIFFERENTIAL MEDIUM FOR *CANDIDA* SPP. 1925

TABLE 1. Growth and colony colors of 726 yeast isolates incubated for 2 days on CHROMagar Candida at 37°C

Species	Total no. of isolates	No. of isolates <3 mo old	Range of colony colors described	Color reference number(s) ^a
<i>Candida albicans</i>	285	117	Green	3258/338
<i>Candida famata</i>	15	12	White, light pink, pink	435/436
<i>Candida (Torulopsis) glabrata</i>	84	44	White, pink, purple	5135/5155
<i>Candida guilliermondii</i>	13	2	Pale pink, purple	435/436
<i>Candida humicola</i>	3	0	Light gray, pink	3/435
<i>Candida inconspicua</i>	3	0	White, pink	435
<i>Candida kefyr</i>	15	0	Pink, purple	435/257
<i>Candida krusei</i>	43	27	Pale pink, purple (rough with spreading, pale edges)	435/5025
<i>Candida lambica</i>	5	0	Pink	5135
<i>Candida lusitanae</i>	21	4	Pink, grayish purple	5135/5155
<i>Candida norvegensis</i>	5	0	White, pale pink (slightly rough with spreading, pale edges)	435
<i>Candida parapsilosis</i>	74	2	White, pale pink	435
<i>Candida pelliculosa</i>	7	0	Pink, pale purple	434
<i>Candida tropicalis</i>	54	6	Dark blue to blue-gray, with dark halo in agar	548/549
<i>Candida utilis</i>	2	0	Pink, purple	434/5155
<i>Cryptococcus laurentii</i>	1	0	Pink-purple	434
<i>Cryptococcus neoformans</i>	23	1	Gray, pale pink	2/434
<i>Debaryomyces polymorphus</i>	2	0	Pink	434
<i>Geotrichum</i> spp.	5	1	Pale to pink (green in agar in two cases) or purple, fimbriate	5405
<i>Pichia</i> spp.	10	1	Highly variable: pink, gray-purple, purple with green edges, dark agar halo in two strains	
<i>Saccharomyces cerevisiae</i>	22	6	White-purple	435/5155
<i>Trichosporon</i> spp.	34	10	Variable; small, pale, "dirty pink" or "dirty gray-green" (becoming darker and rough on prolonged incubation)	5165/5635-5645

^a Color reference numbers are from the Pantone color guide. Colonies of yeasts that gave colonies described as "pale to pink to purple" usually had a dirty, brownish color and a narrow, paler edge (Fig. 1).

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● CHROMagar™
Orientation

JOURNAL OF CLINICAL MICROBIOLOGY, July 1996, p. 1788-1793
 0095-1137/96/\$04.00+0
 Copyright © 1996, American Society for Microbiology

Vol. 34, No. 7

Evaluation of CHROMagar Orientation for Differentiation and Presumptive Identification of Gram-Negative Bacilli and *Enterococcus* Species

JOHN MERLINO,* STEVEN SIARAKAS, GRAHAM J. ROBERTSON, GLENN R. FUNNELL, THOMAS GOTTLIEB, AND ROSS BRADBURY

Department of Microbiology and Infectious Diseases, Concord Repatriation General Hospital, Concord, New South Wales 2139, Australia

Received 14 September 1995/Returned for modification 9 November 1995/Accepted 8 April 1996

Vol. 34, 1996

EVALUATION OF CHROMOGENIC MEDIUM 1791


FIG. 1. Representative chromogenic reactions of selected organisms on CHROMagar. 1, *E. coli*; 2, *P. vulgaris*; 3, *P. mirabilis*; 4, *E. cloacae*; 5, *C. diversus*; 6, *Klebsiella pneumoniae*; 7, *Subdoligranulum typhimurium*; 8, *Aeromonas hydrophila*; 9, *Shigella sonnei*; 10, *P. aeruginosa*; 11, *Y. enterocolitica*; 12, *M. morganii*; 13, *Plesiomonas shigelloides*; 14, *Serratia marcescens*; 15, *A. calcoaceticus*; 16, *Stenotrophomonas maltophilia*; 17, *Enterococcus faecalis*.


FIG. 2. Colonies plated from a mixed suspension, showing two different gram-negative isolates on CHROMagar. The two species can be distinguished by their color appearance. Pink-red, *E. coli*; blue, *Klebsiella pneumoniae*.

tests increased the overall accuracy of identification from 60 to 98% (Tables 1 and 2). Results compared favorably with those obtained by a proven biochemical identification method (11). In comparison with other *E. coli* strains tested, 99.3% gave good correlation with the chromogenic agar. The remaining 0.7% were ONPG-negative strains and were not differentiated by this medium and remained colorless. However, these iso-

1790 MERLINO ET AL.

J. CLIN. MICROB.



TABLE 1. Pigment reaction and presumptive identification of gram-negative and -positive bacteria on CHROMagar Orientation

Organism	Total no. tested	No. (%) with described color	Description of pigment	Pattern color when colonies are*	No. of isolates		
					With correct presumptive identification†	Correctly identified with additional test‡	Not identified by CHROMagar Orientation§
<i>Aeromonas hydrophila</i>	16	16	White nontransparent; initially clear to pink at 37°C (changing to blue with clear edges at RT)	21C	26		16
<i>Alcaligenes denitrificans</i>	1	1	Clear	20C			1
<i>Campylobacter jejuni</i>	7	7	Light bluish green with ill-defined edges	21C		7	
<i>Campylobacter jejuni</i>	15	10 (67)	Violet purple colonies with pinkish border	20C/20C/20C		15	
<i>Enterobacter aerogenes</i>	15	15	Metallic blue with strong purple pink halo (yellowed base)	20C		15	
<i>Enterobacter agglomerans</i>	2	2	Metallic blue	20C		2	
<i>Enterobacter cloacae</i>	32	41 (95)	Metallic bluish purple with strong purple to pink halo around periphery	20C			
<i>Enterobacter cloacae</i>		4 (12)	Metallic bluish purple with no halo	20C			
<i>Enterobacter isolatus</i>	1	1	Metallic blue with strong purple halo	20C/20C		1	
<i>Enterococcus faecalis</i>	46	46	Diffuse blue pinkish within agar (-24 h)	30C			46
<i>Enterococcus faecium</i>	4	4	Diffuse blue pinkish within agar (-24 h)	30C			4
<i>Enterococcus spp.</i>	18	18	Diffuse blue pinkish within agar (-24 h)	30C	18		
<i>Escherichia coli</i>	90	90 (99.3)	Pinkish to red	18C/18C	90		
<i>Escherichia coli</i>		4 (0.7)	Clear	20C		4	
<i>Klebsiella pneumoniae</i>	18	18	Metallic blue with no pink halo	20C		18	
<i>Klebsiella pneumoniae</i>	28	169 (74)	Metallic blue with slight pink halo around periphery	20C		28	
<i>Klebsiella pneumoniae</i>		39 (26)	Metallic blue with no pink halo	20C			
<i>Morganella morganii</i>	26	26	Clear transparent—slight brownish background	30C/30C		26	
<i>Plesiomonas shigelloides</i>	1	1	Transparent pink	30C			1
<i>Proteus mirabilis</i>	14	14	Clear with diffusible brown	30C		14	
<i>Proteus vulgaris</i>	16	11 (69)	Brown green	30C	15		
<i>Proteus vulgaris</i>		1 (6)	Clear with diffusible brown	30C		1	
<i>Providencia stuartii</i>	11	11	Transparent, yellow to green-tinted edges, ill-defined	30C	11		
<i>Providencia stuartii</i>	2	2	Transparent, yellow	Not referenced			2
Group B	2	2	Clear	Not referenced			2
<i>paratyphi</i>	1	1	Clear	Not referenced			1
<i>paratyphi</i>	1	1	Clear	Not referenced			1
<i>typhimurium</i>	1	1	Clear	Not referenced			1
<i>typhimurium</i>	4	4	Light agar blue—may darken to navy blue with pigment production at RT	20C	4		
<i>Serratia marcescens</i>	29	29		20C/30C	29		
<i>Shigella boydii</i>	1	1	Clear pale grey	Not referenced			1
<i>Shigella sonnei</i>	1	1	Pink	18C/18C			1
<i>Stenotrophomonas maltophilia</i>	5	3	Clear	Not referenced			3
<i>Yersinia enterocolitica</i>	1	1	Diffused blue at 37°C	30C			1


* Pattern Color Guide (17) (1988-1990).
 † Presumptive identification based on color and morphology on CHROMagar Orientation only.
 ‡ Identification with additional tests or direct rapid tests as described in Table 2.
 § Could be identified only with a commercial kit or a larger number of subtests.
 ¶ Clear and transparent colonies that could not be referenced by the color guide (RT, room temperature (-24°C)).

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National multidrug-resistant bacteria (MDRB) surveillance: 8-year experience.

Carbonne A, Arnaud J, Maugat S, Marty N, Dumartin C, Bertrand X, Bajolet O, Sarrailh A, Astagneau P, Jarlier V; MDRB Surveillance National Steering Group (BMR-Rai). CCLIN Paris Nord, rue Didot, 75014 Paris, France. anne.carbonne@egp.aphp.fr

Abstract

BACKGROUND: In the mid-1990s, the prevalence rate of multidrug-resistant bacteria (MDRB) spread then became a major priority in the national infection control programme (ICP).

METHODS: To evaluate the impact of the ICP, a national coordination of MDRB surveillance was set up in 2002. Data were collected 3 months a year in healthcare facilities (HCFs) on a voluntary basis. All clinical specimens of methicillin-resistant *Staphylococcus aureus* (MRSA) and extended-spectrum β -lactamase-producing Enterobacteriaceae (ESBLE) were prospectively included. Incidences per 1000 patient days (PDs) were calculated and trends in incidence from 2003 to 2010 were assessed.

RESULTS: Participation in the surveillance increased from 478 HCFs in 2002 to 933 in 2010. In 2010, MRSA incidence was 0.40/1000 PDs: 1.14 in intensive care units (ICUs), 0.48 in acute care facilities (ACFs) and 0.27 in rehabilitation and long-term care facilities (RLTCFs). ESBLE incidence was 0.39/1000 PDs: 1.63 in ICUs, 0.46 in ACFs and 0.23 in RLTCFs. MRSA incidence significantly decreased from 0.72/1000 PDs in 2003 to 0.41/1000 PDs in 2010 ($P < 10^{-3}$); in contrast, ESBLE incidence significantly increased from 0.17/1000 PDs to 0.48/1000 PDs ($P < 10^{-3}$). The most prevalent ESBLE were *Enterobacter aerogenes* (34%) and *Escherichia coli* (25%) in 2003 and *E. coli* (60%) and *Klebsiella pneumoniae* (18%) in 2010.

CONCLUSION: These results demonstrate the positive impact of the national ICP on MRSA rates. In contrast, ESBLE incidence, especially ESBL-producing *E. coli*, is increasing dramatically and represents a serious threat for hospitals and for the community that deserves specific control actions.

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

CDC. Laboratory practices for prenatal group B streptococcal screening and reporting. Connecticut, Georgia, and Minnesota, 1997-1998. MMWR 1999, 48, 426-428.

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MMWR

Morbidity and Mortality Weekly Report

www.cdc.gov/mmwr

Recommendations and Reports

November 19, 2010 / Vol. 59 / No. RR-10

Prevention of Perinatal Group B Streptococcal Disease

Revised Guidelines from CDC, 2010


Continuing Education Examination available at <http://www.cdc.gov/mmwr/cme/conted.html>

DEPARTMENT OF HEALTH AND HUMAN SERVICES
 CENTERS FOR DISEASE CONTROL AND PREVENTION

- Pregnant women should undergo vaginal-rectal screening for GBS colonization at 35-37 weeks.
- Intrapartum antibiotic prophylaxis (IAP) is recommended for:
 - Women who delivered a previous infant with GBS disease
 - Women with GBS bacteriuria in the current pregnancy
 - Women with a GBS-positive screening result in the current pregnancy
 - Women with unknown GBS status who deliver at less than 37 weeks' gestation, have an intrapartum temperature of 100.4°F or greater, or have rupture of membranes for 18 hours or longer.
- Penicillin remains the preferred agent with ampicillin an acceptable alternative.

FIGURE 1. Incidence of early- and late-onset invasive group B streptococcal (GBS) disease — Active Bacterial Core surveillance areas, 1990–2008, and activities for prevention of GBS disease

Abbreviations: ACOG = American College of Obstetricians and Gynecologists and AAP = American Academy of Pediatrics.

Source: Adapted from Jordan HT, Farley MM, Craig A, et al. Revisiting the need for vaccine prevention of late-onset neonatal group B streptococcal disease. *Pediatr Infect Dis J* 2008;27:1057–64.

* Incidence rates for 2008 are preliminary because the live birth denominator has not been finalized.

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Home > Food Safety > Fact sheet > Campylobacter

Main Menu

- About the Agency
- Infectious Diseases
- Chronic Diseases
- Travel Health
- Food Safety
- Immunization & Vaccines
- Emergency Preparedness & Response
- Health Promotion
- Injury Prevention
- Lab Biosafety & Biosecurity
- Surveillance
- Explore
 - Media Room
 - Acts & Regulations
 - Reports & Publications
 - A-Z Index
- Transparency
 - External Advisory Bodies
 - Completed Access to Information Requests
 - Proactive Disclosure

+/- TEXT PRINT SHARE

Campylobacter

Campylobacter refers to a group of bacteria that is commonly found in the intestines of animals, including poultry, wild birds, swine, cattle, rodents and household pets. *Campylobacter jejuni* and *Campylobacter coli* are the two species that most often infect humans, causing an illness called campylobacteriosis.

A person becomes infected from eating or drinking contaminated food, water, unpasteurized or raw milk or from close contact with infected animals or infants.

Campylobacteriosis attacks the digestive system, causing diarrhea, abdominal pain, malaise, fever, nausea and vomiting.

Use Left / Right Arrows to move through menu

Causes Symptoms Risks Treatment Prevention Surveillance

Surveillance

Does the Public Health Agency of Canada keep track of cases of campylobacteriosis across the country?

Yes. The Agency works with the provinces and territories to track the number of laboratory-confirmed cases of campylobacteriosis across the country.

Are campylobacteriosis outbreaks common in Canada?

Campylobacteriosis is the most frequently reported food-related illness in Canada. Many of these illnesses are sporadic cases, but some are part of outbreaks.

Outbreaks of campylobacteriosis are uncommon in Canada, but do occur.

A case refers to illness in one person whereas an outbreak refers to two or more people linked by a common exposure within a specific time frame.

A national outbreak occurs when illness is linked in two or more provinces or territories.

CHROMagar™
Campylobacter

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CHROMagar™
Campylobacter

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DRUG-RESISTANT CAMPYLOBACTER

THREAT LEVEL
SERIOUS


This bacteria is a serious concern and requires prompt and sustained action to ensure the problem does not grow.

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Copyright © 1994, American Society for Microbiology